

LASER PEN

LA-X P500

THE EXPERT'S LASER
FOR THE MOST
EXACTING DEMANDS



THE LASER PEN SETS NEW STANDARDS



The Laser Pen LA-X P500 is suitable for small areas and targeted punctual laser therapy – especially for the treatment of pain, trigger and acupuncture points.

- Point-precise and powerful
- Optimal stimulatory effect and penetration depth
- CW, multi-frequency and frequency series pre-programmed
 - 13 FI bands and 5 sweeps in autorun
 - Free choice of frequencies with additional memory locations
 - Simple user guidance and optimum treatment comfort
 - Laser parameters can be set
 - With point search function and automatic laser
 - Replaceable fibre-optic cables for special applications

LASER PEN

LA-X P500

- 808 nm – 500 mW
10 laser output settings (10 % – 100 %)
- Types of operation: CW, multi-frequency, pre-programmed frequency series (f. i. Nogier, Bahr, Reininger, Fi-Bands [Autorun], Sweeps)
- Freely programmable frequencies
1 to 9.999 Hz / 20 memory locations
- Therapy duration or dose can be set
- Automatic calculation of the therapy time when the dose has been pre-selected
- Memory function of the values saved last (output, frequency))
- Special head for replaceable fibre-optic cables + 1 fibre-optic cable
- Weight 200 g – dimensions 210 x 35 x 32 mm



After the successful market launch of this application method more than 2 decades ago, the LA-X P500 has now become the ultimate technological reference.

Many additional functions, a large colour display, improved functions for routine procedures as well as its great user friendliness are convincing arguments in favour of the new laser system.

From offering the choice between ten laser power settings, to predefined frequency programmes, Fi-Bands, sweeps and freely selectable frequencies,

all the way to conducting high-definition searches for acupuncture points, this unit literally hands the therapist an optimum spectrum of modern laser technology.



INTERESTED?

Request information material at www.mkw-laser.de

MKW Lasersystem GmbH | Landstraße 67 | D-76547 Sinzheim/Baden-Baden
Tel: 0 72 21 - 9 88 391-2 | Fax: 0 72 21 - 9 88 393 | mkw@mkw-laser.de | www.mkw-laser.de

MKW
LASERSYSTEM